



ORTHOPAEDICS INC.

Design and Manufacture of Ponseti Clubfoot Brace

P.O. Box 125
604 North Parkway Street
Wayland, IA 52654
1-877-766-7384

Measurement and size table for the Ponseti AFO. (All measurements are in centimeters.)

Length	Width	Size
7.3	3.8	0000
7.8	4.0	000N
7.8	4.5	000
7.8	5.0	000W
8.3	4.0	00N
8.3	4.5	00
8.3	5.0	00W
9.0	4.0	0N
9.0	5.0	0
9.0	5.5	0W
9.5	4.5	1N
9.5	5.0	1
9.5	5.5	1W
10.5	4.7	2N
10.5	5.2	2
10.5	6.0	2W
11.5	5.5	3N
11.5	6.0	3
11.5	6.5	3W
12.5	6.0	4N
12.5	6.5	4
13.2	5.5	5N
13.2	6.5	5
14.5	6.0	6N
14.5	6.5	6
15.5	6.5	7N
15.5	7.5	7
16.5	7.7	8
17.5	8.0	9
18.5	8.4	10
19.5	8.7	11
20.5	9.0	12

* Please be sure to order the next larger size if measurements fall in-between sizes.
Please contact us if special conditions are present or if measurements do not fall in a provided size range.